

A novel LC/MS/MS method with dual ion source for analysis of 102 pesticides and 5 mycotoxins in hemp



Avinash Dalmia
Senior Principal Application Scientist
PerkinElmer

Outline

- **Pesticides and their action limits in USA states and Canada for hemp/cannabis**
- **Analysis of non-polar pesticides using novel LC/MS/MS method with dual ESI and APCI source**
- **Mitigation of matrix effects with LC and MS method optimization to improve method sensitivity**
- **Sample preparation method – simple solvent extraction**
- **Internal standards – method accuracy**

Oregon, California, Canada & Colorado Action Limits for Pesticides in Hemp/Cannabis

Analyte	OR/ppm	CA/ppm	Canada/ppm	CO/ppm
Acequinocyl	2	0.1	0.03	0.03
Chlorfenapyr	1	0.1	0.05	0.05
Cyhalothrin	NA	NA	NA	0.25
Endosulfan-beta	NA	NA	0.05	0.05
Methoprene	NA	NA	2	2
Quintozene	NA	0.1	0.02	0.02

LC/MS/MS System with Dual ESI and APCI Source


Electrospray Ion (ESI)

Atmospheric Pressure Chemical Ionization (APCI)

LC/MS/MS System




Dual ESI & APCI Source




Non-Polar Pesticides Analyzed with ESI Source



Cyhalothrin



Methoprene



Cyhalothrin in Hemp using LC/MS/MS with ESI Source



Cyhalothrin – 0.1 ppm

S/N = 28


LOQ = 0.05 ppm

Method LOQ of 0.05 ppm meets action limits of 0.25 ppm in hemp

**THE EMERALD
CONFERENCE**

Produced by MJBizScience


Methoprene in Hemp using LC/MS/MS with ESI Source




Method LOQ of 0.1 ppm meets action limits of 2 ppm in hemp

Non-Polar Pesticides Analyzed with APCI Source



Etridiazole



Quintozene (PCNB)



Etridiazole in Hemp with APCI Source



Etridiazole – 0.01 ppm

S/N = 18

LOQ = 0.01 ppm


Method LOQ of 0.01 ppm meets action limits of 0.03 ppm in Hemp

**THE EMERALD
CONFERENCE**

Produced by MJBizScience

© 2023 MJBiz, a division of Emerald X, LLC


Quintozene in Hemp using LC/MS/MS with APCI Source



Quintozene – 0.01 ppm



S/N = 23

LOQ = 0.005 ppm




Method LOQ of 0.005 ppm meets action limits of 0.02 ppm in Hemp

Ionization Mechanism with APCI Source in Negative Mode



Hemp Matrix is Challenging for Pesticides Analysis

- **Complex matrix**
- **Broad Interference from cannabinoids present in 10-20% (100,000-200,000 ppm) range**
- **Terpenes and other compounds are present in high ppm range (10-5000 ppm)**





Propiconazole (0.02 ppm) in Hemp with ESI Source

Matrix Interference




$[M+H]^+ = 342 \rightarrow 69$



Minimal Matrix Interference



$[M+H+2]^+ = 344 \rightarrow 69$




Green trace is signal for 20 ppb propiconazole in hemp


Red trace is signal for hemp blank

Endosulfan- β (0.1 ppm) in Hemp with APCI Source

Fast 6 minute LC Gradient
75 % Ion Suppression



Optimum 11 minute LC Gradient
Insignificant Ion Suppression



Method LOQ of 0.01 ppm meets Action Limits of 0.05 ppm in Hemp


**THE EMERALD
CONFERENCE**

Produced by MJBizScience


© 2023 MJBiz, a division of Emerald X, LLC

MS Scan for CBDA Retention time with 2 LC Methods


Fast 6 min LC Gradient
Endosulfan- β , RT = 3.28 min
CBDA, RT = 3.29 min




Optimum 11 min LC Gradient
Endosulfan- β , RT = 3.93 min
CBDA, RT = 4.03 min




Acequinocyl in Hemp with ESI & APCI Source



ESI Source (0.01 ppm)



APCI Source (0.01 ppm)



Method LOQ of 0.005 ppm meets Action Limits of 0.03 ppm in Hemp

Sample Preparation Method - Solvent Extraction

0.5 gm of Ground Hemp in 50 ml tube



2.5 ml of ACN with 0.1 % formic acid (FA)



Vortex 10 min



Centrifuge 5 min at 3000 rpm




Dilute 500 μ L extract with 490 μ L of ACN/0.1 % FA and 10 μ L of IS

Sample Preparation Methods

- Solvent Extraction – (ACN extraction)
- Solvent Extraction with Evaporation
- QuEChERS
- Dispersive Solid Phase Extraction (dSPE)
- Solid Phase Extraction (SPE)
- Immunoaffinity Columns

Simple, Speed, Cheap, Green,
Automation & Extraction
Efficiency



Automated Testing Workflow

FULLY AUTOMATED CANNABIS & HEMP PESTICIDE-RESIDUE TESTING WORKFLOW



WEIGHT & SAMPLE ID

SAMPLE ID & WEIGHT


SAMPLE LIST & BATCH
INFORMATION

BATCH & QUALITY
CONTROL INFORMATION

SAMPLE RESULTS &
BATCH INFORMATION

SAMPLE & QC
TEST DATA

METRC &
COA Report



RECEIVE & WEIGH SAMPLES

Measure out representative Cannabis or Hemp sample




SAMPLE PREPARATION

JANUS[®] G3 420 Workstation with ONE Pesticide^{420™} reagent kit, vortex, centrifuge & integrated barcode scanner



INSTRUMENT ANALYSIS

QSign[®] 420 LC/MSMS System - Completes the entire pesticide method



DATA PROCESSING

Full Regulatory Compliance & Client COA Report


THE EMERALD
CONFERENCE

Produced by MJBizScience

© 2023 MJBiz, a division of Emerald X, LLC

$$\text{Extraction Efficiency} = \frac{\text{Prespiked extract} * 100}{\text{Postspiked Extract}}$$

Extraction Efficiency with Solvent Extraction



Extraction Efficiency – Different Sample Preparation

Analyte	QuEChERS with dSPE(1)	Solvent Extraction/SPE/ dSPE(2)	Solvent Extraction/ Evaporation(3)	Solvent Extraction
Abamectin	ND	18%	ND	94%
Fenpyroximate	74%	61%	44%	97%
Imazalil	49%	21%	114%	93%
Pyridaben	68%	58%	ND	95%
Spirotetramat	76%	63%	132%	98%


1. www.sigmaaldrich.com

2. Moulins, et al., J. AOAC International, (2021) vol. 101, 1948-1960.

3. Michlig, et al., J. Chrom. A (2021), doi:<https://doi.org/10.1016/j.chroma.2021.462097>

$$\text{Matrix Effect} = \frac{(\text{Postspiked extract} - \text{solvent standard}) * 100}{\text{Solvent standard}}$$

Hemp Matrix Effect



Internal Standards – Accounting Matrix Effects

Internal Standard

Acequinocyl-d25

Atrazine-d5

Bifenthrin- d5

Boscalid-d4

Captan-d6

Carbaryl-d7

Chlorpyrifos (Diethyl-d10)

Coumaphos-d10

Daminozide-d4

Diazinon-d10

Dichlorvos-d6

Dimethoate-d6

Doramectin

Fenoxycarb-d3

Fipronil-¹³C2, ¹⁵N2

Internal Standard

Fludioxonil-¹³C3

Imidacloprid-d4

Malathion-d6

Methylparathion-d6

Myclobutanil-d9

Naled-d6

Ochratoxin-A-d5

Pentachloronitrobenzene-¹²C6

Permethrin-d5

Phosmet-d6

Piperonylbutoxide-d9

Pyridaben-d13


Spirotetramat-d5

Thiamethoxam-d3

Trifloxystrobin-d6

$$\text{Accuracy} = \frac{\text{Measured Concentration} * 100}{\text{Spiked Concentration}}$$

Method Accuracy (Internal Standard)



Method LOQ and AL for Pesticides & Mycotoxins in Hemp

Analyte	LOQ(AL)/ppm	Analyte	LOQ(AL)/ppm	Analyte	LOQ(AL)/ppm	Analyte	LOQ(AL)/ppm
Abamectin	0.02(0.1)	Daminozide	0.01(0.1)	Hexythiazox*	0.005(0.01)	Propiconazole	0.02(0.1)
Acephate	0.005(0.02)	Deltamethrin	0.02(0.5)	Imazalil	0.01(0.05)	Propoxur	0.005(0.02)
Acequinocyl*	0.005(0.03)	Diazinon	0.005(0.02)	Imidacloprid	0.005(0.02)	Pyraclostrobin	0.005(0.02)
Acetamiprid	0.005(0.1)	Dichlorvos	0.01(0.1)	Iprodione*	0.01(1)	Pyrethrins	0.02(0.05)
Aldicarb	0.01(1)	Dimethoate	0.005(0.02)	Kinoprene	0.1(0.5)	Pyridaben	0.005(0.05)
Allethrin	0.05(0.2)	Dimethomorph	0.01(0.05)	Kresoxim-methyl	0.01(0.02)	Pyriproxyfen	0.005(0.01)
Atrazine	0.005(0.025)	Dinotefuran	0.005(0.1)	(Lambda) Cyhalothrin	0.05(0.25)	Quintozene*	0.01(0.02)
Azadirachtin	0.05(1)	Diuron	0.005(0.125)	Malathion	0.01(0.02)	Resmethrin	0.02(0.1)
Azoxystrobin	0.005(0.02)	Dodemorph	0.005(0.05)	Metalaxyl	0.005(0.02)	Spimetroram	0.01(0.02)
Benzovindiflupyr	0.005(0.02)	Endosulfan sulfate*	0.01(0.05)	Methiocarb	0.01(0.02)	Spinosad	0.01(0.1)
Bifenazate	0.005(0.02)	Endosulfan-alpha*	0.1(0.2)	Methomyl	0.005(0.05)	Spirodiclofen	0.02(0.25)
Bifenthrin	0.01(1)	Endosulfan-beta*	0.01(0.05)	Methoprene	0.2(2)	Spiromesifen	0.01(3)
Boscalid	0.01(0.02)	Ethoprophos	0.005(0.02)	Mevinphos	0.01(0.05)	Spirotetramat	0.01(0.02)
Buprofezin	0.005(0.02)	Etofenprox	0.01(0.05)	MGK-264	0.01(0.05)	Spiroxamine	0.005(0.1)
Carbaryl	0.01(0.05)	Etoxazole	0.005(0.02)	Myclobutanil	0.005(0.02)	Tebuconazole	0.005(0.05)
Carbofuran	0.005(0.02)	Etridiazole*	0.01(0.03)	Naled	0.01(0.1)	Tebufenozide	0.005(0.02)
Chlorantraniliprole	0.005(0.02)	Fenhexamid	0.01(0.125)	Novaluron	0.005(0.05)	Teflubenzuron	0.01(0.05)
Chlorfenapyr*	0.02(0.05)	Fenoxycarb	0.005(0.02)	Oxamyl	0.005(3)	Tetrachlorvinphos	0.005(0.02)
Chlorpyrifos*	0.01(0.04)	Fenpyroximate	0.005(0.02)	Paclobutrazol	0.01(0.02)	Tetramethrin	0.02(0.1)
Clofentezine	0.005(0.02)	Fensulfthion	0.01(0.02)	Parathion-methyl*	0.005(0.05)	Thiabendazole	0.01(0.02)
Clothianidin	0.005(0.05)	Fenthion	0.005(0.02)	Permethrin	0.01(0.5)	Thiacloprid	0.005(0.02)
Coumaphos	0.005(0.02)	Fenvalerate*	0.02(0.1)	Phenothrin	0.02(0.05)	Thiamethoxam	0.005(0.02)
Cyantraniliprole	0.005(0.02)	Fipronil*	0.005(0.06)	Phosmet	0.01(0.02)	Thiophanate-methyl	0.01(0.05)
Cyfluthrin	0.2(0.2)	Fonicamid	0.005(0.05)	Piperonyl butoxide	0.02(0.2)	Trifloxystrobin	0.005(0.02)
Cypermethrin	0.1(0.3)	Fludioxonil*	0.005(0.02)	Pirimicarb	0.005(0.02)	Aflatoxin(B1+B2+G1+G2)	0.012(0.020)
Cyprodonil	0.01(0.25)	Fluopyram	0.005(0.02)	Prallethrin	0.01(0.05)	Ochratoxin A	0.004(0.005)

Conclusions

- **LC-MS/MS method with dual APCI and ESI source analyzes all of pesticides and mycotoxins regulated by different USA states and Canada to meet their action limits in hemp**
- **Matrix interference and matrix effects were minimized by MS and LC method optimization and by changing mode of ionization to improve method sensitivity**
- **Solvent extraction can achieve high throughput with automation, less cost and excellent extraction efficiency**
- **Internal standards compensated for matrix effects such as ion suppression/enhancement to get overall accuracy of 70-120 %**

Acknowledgements

- **Jacob Jalali, Erasmus Cudjoe, Saba Hariri, Sheng-Suan Cai, Feng Qin, Hui Qiao, Alex Kasperkiewicz, Scott LaPaglia, Jeff Wu, Tyrally Ordinario, Toby Astill, Charlie Schmidt, Tony Rhoden, John Rontree (PerkinElmer)**
- **Charles Johnson, Joey Kingstad, Jeremy Campbell, Kevin Smith, Ben Armstrong, Luke Ward, Adam Floyd, Molly Murphy, Paul Dorenbach & Travis Ruthenburg (Cannabis testing Labs)**

Thank you for your Attention!

Avinash Dalmia

Senior Principal Application Scientist, PerkinElmer

Avinash.Dalmia@perkinelmer.com

203-558-1748

<http://www.perkinelmer.com/cannabis>

**THE EMERALD
CONFERENCE**

Produced by MJBizScience

© 2023 MJBiz, a division of Emerald X, LLC

We want to hear from you!

Scan the QR code below to provide your feedback on the presentation.

