

Comprehensive comparison of cannabis cultivars using untargeted analytical workflows

Aldwin Anterola¹, Laura McGregor² and John Abrams³


¹Chief Scientific Officer, Veda Scientific

²SepSolve Analytical

³CESC

Cannabis aroma volatiles

- Cannabis terpenes and terpenoids
 - ✓ Impart distinctive aromas and flavours




"From the literature, it is therefore difficult to know what to believe when it comes to the odor of limonene." (J. Chem. Educ. 2021, 98, 11, 3600–3607)

Many other aroma volatiles are also present, e.g. esters, aldehydes, etc.

Aroma volatiles vs. cultivars


 Much remains unknown about cultivar distinctions

 Specific aroma profiles can be engineered by plant breeders




Source: <u>leafly.com</u>

Is GC-MS fit-for-purpose?




Is GC-MS fit-for-purpose?






Challenges in terpenes analysis

- Over 200 terpenes have been identified in cannabis
- 1D GC results in coelutions between terpenes and terpenoids
- Abundance of important terpenes can be over-estimated





What is GC×GC?





What is GC×GC?






Cannabis terpenes by GC×GC-TOF MS

 Enhanced separation of terpenes and terpenoids

 Comprehensive detail on terpene profile

Detection of other compounds






Analysis of cannabis cultivars

- Dried cannabis flowers were obtained from 12 different cultivars
- Extracted with methanol
- Analyzed beta-pine/ limonene ratio
- Looked at other features


Cultivar	BP/LIM Ratio
Pineaple Upside Down Cake	Hi
Blue Dreams	Hi
Northern Wreck	Lo
Acapulco Gold	Hi
P-Funk	Hi
P-Funk x Black Dog	Hi
Lemon OG	Lo
Zkittlez x Black Dog	Hi
Z-Dog	Lo
Orange Valley Delight	Hi
Super Sour Diesel	Lo
GMO x Black Dog	Lo

Diagnostic ratios in cannabis profiling






Diagnostic ratios in cannabis profiling




Diagnostic ratios in cannabis profiling



β-Caryophyllene/α-HumuleneRatio


Analysis of cannabis by GC×GC-TOF MS



More detail than a single ratio can provide




Reducing the dimensionality of data



• Untargeted, tile-based data analysis

Expanding on the Floral vs. Fuel classification:




Categorization based on beta-pinene / limonene ratio type


- 12 cannabis cultivars analysed in triplicate
- Untargeted, tile-based data analysis




Ratio of 2 BCP Features Correlates Largely with BPIN / LIM Ratio Type



Expanding on the Floral vs. Fuel classification: Categorization based on beta-pinene / limonene ratio type




Distinguishing individual cultivars




- All cultivar replicates cluster closely
- Z-Dog (ZD) is clearly differentiated by PC1
- Super Sour Diesel (SSD) is clearly differentiated by PC2





Distinguishing individual cultivars






Distinguishing individual cultivars








Looking beyond terpenes




• Feature identified as the ester, ethyl 3-methyl-2-butenoate




Looking beyond terpenes






¹Good Scents database, http://www.thegoodscentscompany.com/


α -Pinene primarily appears as a cluster of features





α -Pinene feature cluster includes sabinene . . .







β-Ocimene-15.940_0.375_94-15.94(0.00625) (18) β-Ocimene-15.940_1.575_93-15.94(0.02625) (18)


Summary


- GC×GC—TOF MS provides the enhanced separation necessary for robust profiling of terpenes
- Additional chemical differences were uncovered between the "Floral" and "Fuel" classes
- Unique chemical markers can further distinguish individual cultivars
- Provides improved product definition, through superior cultivar selection and improved labelling and categorization strategies



Cannabis aroma volatiles

- Cannabis terpenes and terpenoids
 - ✓ Impart distinctive aromas and flavours





"From the literature, it is therefore difficult to know what to believe when it comes to the odor of limonene." (J. Chem. Educ. 2021, 98, 11, 3600–3607)

Many other aroma volatiles are also present, e.g. esters, aldehydes, etc.

Thanks for listening! Any questions?



aldwin@vedascientific.co



+1 805-741-8881



We want to hear from you!

Scan the QR code below to provide your feedback on the presentation.

